




AEROSOL
MAGEE SCIENTIFIC

SIMPLE AND STRAIGHTFORWARD DATA ANALYSIS



CAAT

CARBONACEOUS AEROSOL
ANALYSIS TOOL

Advanced
Black Carbon
Measurement

KEY FEATURES

- IMPORT your data from Magee Scientific instrument
- VALIDATE your data.
- ANALYZE your validated data by using several procedures: Diurnal profile, Time Average, Absorption Coefficient, Ångström exponent, Attenuation, Source Apportionment.
- DISPLAY your data by using various chart options.

APPLICATIONS

- Analyze and display the data from Aerosol Magee Scientific instruments



AEROSOL MAGEE SCIENTIFIC

PRODUCT SPECIFICATIONS

Analyze your Magee Scientific instrument data with CAAT Software.



BENEFITS

Be assured that your data is correct and consistent
Use CAAT as your primary tool for your analysis of data obtained from Magee Scientific instruments.
Immediately check the proper "data status".
Start using CAAT Software tool
No additional requirements.
Straightforward data evaluation and display.
User friendly controls.

EASY INSTALLATION

- Login or Register at [Magee Scientific User's group](#):
Search for CAAT product.
- Go to SOFTWARE tab and download installation files.
- Search for the downloaded files in your computer.
Unzip to a selected destination folder.
- Click on the "CAAT-Setup.exe" and follow the guided procedure.
- You can start using the CAAT application immediately after installation is completed.

CAAT



AEROSOLMAGEESCI.COM

Aerosol d.o.o.
Kamniška 39 A
SI-1000 Ljubljana
Slovenia
+386 1 439 1700

Manufactured in EU by Aerosol d.o.o.

CAAT specification version 1.1 / 02 2023

Specifications are subject
to change without notice.

Keeping an Eye on the Air